





Research services: Biotechnology and chemistry
Research centers: Kobe, Kashiwa, Saitama, Chiba, Shounan iPark (Office)
Leadership:

- Gaku Shimaoka: President and Representative Director, TechnoPro Inc.
- Masami Hayafune: Representative director, TechnoPro R&D Company, TechnoPro Inc.
- Neuroscience and Drug Discovery in public private partnership
  - Revealing neural function in drug discovery with integrated research of neuroscience and electronics



Gene/Cell
Expression vector
transduced cell
line, gene analysis















### **Neuroscience and Drug Discovery**

### Assessment system and applications with;

- Human iPSC-neuron
- MEA: Microelectrode array for detection of extracellular action potential
- AI: Proprietary Artificial intelligence and multivariate analysis for pharmacology

### (1) Human Neurons

- CNS: Cortical neuron, Dopaminergic neuron ...
- Sensory neuron, Motor neuron ...
- Glia cell: Microglia, astrocyte, ...
- **→** Difference of cells from rodent to human

# (2) Efficacy and Pharmacology

- Disease-models: siRNA-KD, gene-mutation for PD, ALS, Alzh, ...
- Patient-derived iPSC neuron: HD, Dravet (Organoid), ...

# (3) Toxicology

- Ion channels/receptors on cortical neuron
- Addiction in dopaminergic neuron
- Adverse event in sensory neuron
- •3D: Application of cerebral organoid



Pharmaceuticals Biotech Academia Research Inst. ...

## **Platform Technology** (Selected)

### **■ CNS**: Seizure, Addiction

- Principal Component Analysis to Distinguish Seizure Liability of Drugs in Human iPSCell-Derived Neurons (Toxicological Sciences, 184:265, 2021)
- Can we panelize seizure? (Toxicological Sciences, 179:3, 2021)
- Toxicological evaluation of convulsant and anticonvulsant drugs in human induced pluripotent stem cellderived cortical neuronal networks using an MEA system (Scientific Reports, 8:10416, 2018)
- Evaluation of drug addiction using human iPSC-dopamine neuron (Japanese Society of Toxicology, 2022)

### **■ Sensory neuron:** Pain assessment

- In vitro pain assay using human iPSC-derived sensory neurons and microelectrode array (Toxicological Sciences, 188:131, 2022)
- Motor neuron: Neuromuscular disease (In preparation for publication)
  - ALS and related diseases models with human iPSC-motor neurons
  - Efficacy for mutated neurons by gene-editing of a disease gene
  - Inhibition effect for accumulation of mutated proteins in patient motor neurons
- AXON tracking: Conduction velocity, synaptic plasticity

  Versatile live-cell activity analysis platform for characterization of neuronal dynamics at single-cell and network level (Nature Communications, 11:4854, 2020)
- **Artificial intelligence**: Raster plots machine learning to predict the seizure liability of drugs and to identify drugs (Scientific Reports, 12:2281, 2022)
- **Drug response**: Physiological maturation and drug responses of human induced pluripotent stem cell-derived cortical neuronal networks in long-term culture (Scientific Reports, 6:26181, 2016)



















# **■** Analysis: Artificial Intelligence



### [Deep learning]

- Image recognition with raster plot of positive controls
- Weighted pattern recognition extracting feature value
- The AI is created to pile up multi-layers of algorithm.


### [Identification]

The AI identifies/classifies feature value between control and test compounds to toxicity risk and efficacy.


### [Advantages]

- A wide range of accurate identification with 4000+ parameters.
- Hi-sensitive evaluation in weak neural firing.

# 1. Separated spike stamp table




### 2. Image Recognition



3. Pattern Recognition


Neural network

#### Pilocarpine, muscarinic R agonist



#### 4. Identification

Kainic acid, glutamate R activator







info-rd@technopro.com
www.technopro.com/en/company/company\_info
(in JPN) www.technopro.com/rd/services/contract/
lineup/pharmacological-evaluation/

Copyright © TechnoPro, Inc. All Rights Reserved.